

# firefly WeGo




Using the GoTo Seat with the WeGo will slightly reduce the seat depth.

**Please note:** The WeGo offers no postural support without the use of the GoTo Seat

To help you get the most out of your WeGo, our clinical team have pulled together some top tips and things that you might find useful.

- 1 If your child has a seat depth measurement of **26cm or below** **A** – adjust the footrest to its highest setting. This will stop your child from sliding forward. **B**



- 2 The GoTo Pressure Relief cushion can be used in conjunction with the GoTo and WeGo together. Folding the foam roll under the cushion will add tilt to the seat, which can help stabilise the pelvis.

**Note: Doing this will reduce the height of the GoTo Seat.**



- 3 Opening up the recline angle, reducing the angle of the footplates or a combination of both will help accommodate children with tight hamstrings.

