




GoTo

Factsheet #06:

**Putting the upholstery back on
and assembling GoTo**

Putting the upholstery back on and assembling GoTo

Seat Base



- 1 Retrieve the straps for the GoTo insert, with straps A, B, C and D at the rear of the GoTo base and straps E and F (the longest straps) at the seat pad.



- 2 Thread the buckle of each black strap (not the fastening) through the holes in the upholstery, followed by the hole in the insert.

Seat Base

- 3** As with the above straps, thread the buckle of the groin strap through its corresponding hole in the upholstery at the bottom of the seat pad, followed by the plastic base.



- 4** Once all of the straps are inserted, place the insert into its upholstery.



Slide the widest part of the GoTo Seat base into the upholstery first. For Size 1, open the Velcro of the seat base before this.

Lateral Support



Before starting, plan which lateral support strap needs to be attached to which side of the upholstery in order for the buckle to fasten.



- 1 Loosen the lateral support straps on the central buckle as much as possible. Align the upholstered cover, plastic lateral support insert and straps as shown above.



- 2 Insert the lateral support strap into the lateral support upholstery.



- 3 Pull the strap through the hole of the upholstery.



Thread buckle through inner hole of upholstery for smaller children.

Lateral Support



- 4 Thread strap through the lateral support insert as shown (this leaves the buckle on the outside of the lateral support to prevent any irritation during use). You can use a different thread pattern to change the width of the lateral to suit your child.




- 5 Place the plastic lateral support insert into the upholstery and pull on the strap to pull the lateral support insert into the upholstery. The central section should appear as shown.



Lateral Support

- 6** Remove the knob that is keeping the slider component in place. Place the lateral support over the slider mechanism component and secure using the knob.



- 7** Turn your GoTo over and clip all the straps together. Check that all the straps are oriented correctly as shown above.
-

Advanced Headrest



1 Overview of parts.




2 Place foam mould onto plastic headrest insert and place halfway into advanced headrest upholstery.


3 Pull the upholstery straps through the foam mould and then through the top holes of the plastic headrest insert (not where the shoulder straps are placed) and secure the poppers.



Advanced Headrest



- 4 Attach shoulder strap buckles to the bottom holes of headrest.



- 5 Align screw hole of headrest to the gap in GoTo base. Fix headrest to GoTo base by screwing knob clockwise.

Standard Headrest



- 1 Insert shoulder strap buckle straight through the front of the headrest insert and secure at the back.



- 2 Slide headrest insert into its upholstery and fasten poppers.



- 3 Attach headrest to GoTo base by aligning screw hole with its gap. Fix to place by screwing the knob clockwise.
-