



**Playpak Factsheet #02:  
How to unpack and  
pack Playpak**

## Unpacking

Playpak arrives ready for use.

Lay it on the floor. Lift up the flap of the bag and unzip the sides by pulling on the orange Firefly zip pulls (watch your fingers in the zips).

Fold flat by gently pulling apart the Velcro on the two short sides.

Remove the shoulder strap using the two hooks on either end before using your Playpak, as it could create an uneven surface underneath.

Now you're ready.

- 1 Rolls (1x long and 1x short)**  
These rolls can be used on their own or with other support elements to provide just the right level of support.
- 2 Positioning Strap**  
Use this versatile strap in addition to the rolls for extra support when needed at the front, sides or back.
- 3 Head/Tummy Support**  
Can be used on its own to support head or trunk, or with other components for extra positioning possibilities.
- 4 Horseshoe Sitting Supports (1x large and 1x small)**  
These profiled cushions are designed to give your child lower and upper back support, while allowing the hands to be free for play and function. When used in conjunction with the other elements, additional front, rear or side support can be increased.



## Packing

There are a number of different ways to pack Playpak efficiently, and you will find the best way for you with practice. This is our recommended technique for efficient packing:

- 1 Use the Velcro on the two small side flaps to stand the flaps up. This will give you an accurate idea of how much space you have.
- 2 Place the large red horseshoe sitting support on to the mat first. The back of the support should be sitting flush with the folding edge of the bag.
- 3 Place the small green horseshoe sitting support in the middle of the large red one.
- 4 The blue head/tummy support will sit neatly on top of the green horseshoe support, along with the positioning strap
- 5 Place the long orange roll at the other folding edge of the bag, as shown, and place the short yellow roll on top of it, (you may need to fold the long orange roll in half)
- 6 Zip up this side of the bag using the orange Firefly zip pull, being careful not to nip your fingers.
- 7 Zip up the other side of the bag, then pull the flap over and push down so that the Velcro 'sticks' securely.
- 8 Attach your shoulder strap using the hooks on either end of the straps.



You are ready to go!