



upsee

Upsee Factsheet #08: Technical Information



Size	Extra Small (1-2 yrs)	Small (2-4 yrs)	Medium (4-6 yrs)	Large (6-8 yrs)
Shoulder to Hip	11"-13.5" (28-32cm)	12"-14" (31-36cm)	14"-15.5" (35-39cm)	15"-16.5" (38-42cm)
Chest Circumference	19"-22" (50-57cm)	20"-24" (52-60cm)	22"-26" (57-65cm)	25"-28" (63-71cm)
Lower Hip Circumference	19"-22" (50-57cm)	20"-24" (52-60cm)	23"-26" (58-67cm)	25"-29" (64-73cm)
Maximum User Weight	33lbs (15kg)	33lbs (15kg)	44lbs (20kg)	55lbs (25kg)

Shoe Size Chart	UK	US	Europe	JP / CN
Infant Sandal	3-7	3.5-7.5	19-24	11.5-15
Standard Child Sandal	5-13	5.5-13	22-32	13.5-20
Adult	3-11	4-12	36-46	23-29.5

Adult Belt (One Size) to fit hip 24" - 42" (60cm - 107cm)

For stability and safety, sandal sizes are automatically allocated to the child harness size selected.

- Infant Sandal is automatically provided with the Extra Small Child Harness
- Standard Child Sandal is automatically provided with Small - Large Child Harness
- Adult Sandal on both options will fit adult shoe 3 - 11 (UK), 4 - 12 (US), 19 - 29 (EU)

Please note that these sizing charts are a guide only. If the child wears AFO's, Splints or body casts, you may need to reconsider which size they require.