

Introducing

splashy & splashy BIG



Messy Play



Out in the garden



At the beach

Splashy & Splashy BIG are portable activity seats which are multi-functional and can help children to participate in a wide range of activities.

These seats include unique components which give children the postural support and comfort that they require. Splashy and Splashy BIG have been designed for outdoor or indoor use, and are lightweight and portable. Splashy has been designed for children aged 1-8 (maximum weight 88 Lbs), whilst Splashy BIG has been designed for children aged 7-14 (maximum weight 121 Lbs). Splashy and Splashy BIG can be used with children on levels III-V of the Gross Motor Function Classification System (GMFCS).

Splashy and Splashy BIG can be used for the following applications

- Messy or sensory play (as floorsitter) at home or school
- Outdoor postural support seat in garden or at beach/park
- As a bath or shower seat

Splashy and Splashy BIG include the following unique components.

Seat cover is warm-touch EVA foam material, which is very comfortable on the skin

Seat can be reclined to 26 different angles to ensure child is in a comfortable position – this is particularly useful for children with limited head control

Support bumpers can be adjusted according to the postural support needs of the child – Eg. Head, lateral, or leg support

Harness can be quickly adjusted from a 3 point (pelvic support) harness to a 5 point harness depending on the individual needs of the child



Splashy BIG

Splashy

Splashy BIG includes a stabiliser to give improved stability when used in the living room, garden or beach

Insurance Guidance

Due to the versatile applications and unique components, Splashy and Splashy BIG should be coded as E1399 for insurance purposes.