



GoTo

The **GoTo Seat** is a safe, supportive, lightweight and portable seat designed to provide additional upper body postural support for children with disabilities. The GoTo comes with a floorsitter and a pressure relief cushion.

Dealer		Marked for		Account	
Date		PO		Ordered by	
Phone		Fax		Email	
Ship to address					

Product Description	Product Codes	Product Chosen	MSRP
GoTo Size 1 - Standard Headrest / Orange / Yellow with floorsitter	167-501-06-SF-CU		\$1,180
GoTo Size 1 - Standard Headrest / Light Blue / Dark Blue with floorsitter	167-501-05-SF-CU		\$1,180
GoTo Size 1 - Standard Headrest / Red Vinyl with floorsitter	167-501-09-SF-CU		\$1,300
GoTo Size 1 - Standard Headrest / Pink / Purple with floorsitter	167-501-04-SF-CU		\$1,180
GoTo Size 1 - Advanced Headrest / Orange / Yellow with floorsitter	167-501-06-CF-CU		\$1,230
GoTo Size 1 - Advanced Headrest / Light Blue / Dark Blue with floorsitter	167-501-05-CF-CU		\$1,230
GoTo Size 1 - Advanced Headrest / Red Vinyl with floorsitter	167-501-09-CF-CU		\$1,350
GoTo Size 1 - Advanced Headrest / Pink / Purple with floorsitter	167-501-04-CF-CU		\$1,230
GoTo Size 2 - Standard Headrest / Orange / Yellow with floorsitter	167-502-06-SF-CU		\$1,380
GoTo Size 2 - Standard Headrest / Light Blue / Dark Blue with floorsitter	167-502-05-SF-CU		\$1,380
GoTo Size 2 - Standard Headrest / Red Vinyl with floorsitter	167-502-09-SF-CU		\$1,500
GoTo Size 2 - Standard Headrest / Pink / Purple with floorsitter	167-502-04-SF-CU		\$1,380
GoTo Size 2 - Advanced Headrest / Orange / Yellow with floorsitter	167-502-06-CF-CU		\$1,450
GoTo Size 2 - Advanced Headrest / Light Blue / Dark Blue with floorsitter	167-502-05-CF-CU		\$1,450
GoTo Size 2 - Advanced Headrest / Red Vinyl with floorsitter	167-502-09-CF-CU		\$1,570
GoTo Size 2 - Advanced Headrest / Pink / Purple with floorsitter	167-502-04-CF-CU		\$1,450
Firefly WeGo - Push Chair Accessory for GoTo Seat	182-500		\$ 1,195
Firefly Backpack - Navy	190-1800		\$178

GoTo Seat	Seat base to top of headrest Max (inches)	Seat base to top of headrest Min (inches)	Lateral support width Max (inches)	Lateral support width Min (inches)	Seat Depth (inches)	Max Weight (lbs)	Age Range (approx)
Size 1	21	14	10	7.5	7	33	1-3 years
Size 2	27	20	12	8	9.5	66	3-8 years