



## GoTo

The **GoTo Seat** is a safe, supportive, lightweight and portable seat designed to provide additional upper body postural support for children with disabilities. The GoTo comes with a floorsitter and a pressure relief cushion.

Dealer		Marked for		Account	
Date		PO		Ordered by	
Phone		Fax		Email	
Ship to address					

Product Description	Product Codes	Product Chosen	MSRP
GoTo Size 1 - Standard Headrest / Light Blue / Dark Blue with floorsitter	167-501-05-SF-CU		\$1,275
GoTo Size 1 - Standard Headrest / Red Vinyl with floorsitter	167-501-09-SF-CU		\$1,405
GoTo Size 1 - Standard Headrest / Pink / Purple with floorsitter	167-501-04-SF-CU		\$1,275
GoTo Size 1 - Advanced Headrest / Light Blue / Dark Blue with floorsitter	167-501-05-CF-CU		\$1,330
GoTo Size 1 - Advanced Headrest / Red Vinyl with floorsitter	167-501-09-CF-CU		\$1,460
GoTo Size 1 - Advanced Headrest / Pink / Purple with floorsitter	167-501-04-CF-CU		\$1,330
GoTo Size 2 - Standard Headrest / Light Blue / Dark Blue with floorsitter	167-502-05-SF-CU		\$1,490
GoTo Size 2 - Standard Headrest / Red Vinyl with floorsitter	167-502-09-SF-CU		\$1,620
GoTo Size 2 - Standard Headrest / Pink / Purple with floorsitter	167-502-04-SF-CU		\$1,490
GoTo Size 2 - Advanced Headrest / Light Blue / Dark Blue with floorsitter	167-502-05-CF-CU		\$1,565
GoTo Size 2 - Advanced Headrest / Red Vinyl with floorsitter	167-502-09-CF-CU		\$1,695
GoTo Size 2 - Advanced Headrest / Pink / Purple with floorsitter	167-502-04-CF-CU		\$1,565
Firefly WeGo - Push Chair Accessory for GoTo Seat	182-500		\$ 1,290
Firefly Backpack - Navy	190-1800		\$195

GoTo Seat	Seat base to top of headrest Max (inches)	Seat base to top of headrest Min (inches)	Lateral support width Max (inches)	Lateral support width Min (inches)	Seat Depth (inches)	Max Weight (lbs)	Age Range (approx)
Size 1	21	14	10	7.5	7	33	1-3 years
Size 2	27	20	12	8	9.5	66	3-8 years